

ERC Catalysts, Inc. 215 Main Street, Miltona, MN 56354 U.S.A. Phone: 1.218.943.7904

LUMBROKINASE (Megadriles earthworm)

SPECIFICATION SHEET CAS # 556743-18-1

DESCRIPTION

Lumbrokinase is a naturally occuring proteoytic and fibrinolytic enzyme derived from the fermentation of Lumbricus rubellus (Megadriles earthworm) from Non-GMO sources. It is best known for supporting cardiovascular health.

PHYSICAL PROPERTIES

A tan to brown free-flowing powder, soluble in water, with characteristic odor and taste.

ENZYMATIC PROPERTIES

The enzyme is stable in a pH range from 6.0 to 10.0 with the optimum pH between 7.0 to 9.0. The optimum temperature is 50°C with a sharp drop over 60°C.

ACTIVITY

One (1) FU unit is defined as the amount of the enzyme which increases the absorbency of the filtrate at 275 nm by 0.01 per minute under the conditions of the assay (JBSL). The acceptance criteria for all enzyme assays is: NLT 85.0% and NMT 115.0 % of the declared units of enzyme activity.*

COUNTRY OF ORIGIN USA

STORAGE/SHELF LIFE/STANDARD PACK SIZE

Product is stable for two years (24 months) if stored at or below 10°C in sealed poly bags in boxes or drums away from sunlight and high humidity. Product is packed in 25 kilo fiber drums or double-wall boxes.

HANDLING PRECAUTIONS

Avoid the formation of aerosol and dust of the product. Repeated inhalation of enzyme aerosol or dust may cause allergic type reactions in sensitized individuals. For detailed information please refer to the SDS.

Description Activity: Identity: Moisture:		Specification NLT 15,000 FU/GM Lumbrokinase NMT 10%	Method JBSL FTIR Ohaus MB-45
Metals: Microbiolo	Lead gical Data:	NMT 0.5 ppm	AOAC 2013.06
	TPC	<10,000 CFU/g	Soleris / AOAC 990.12
	E.coli	Negative/10g	Soleris / AOAC 991.14
	Entero	<100 CFU/g	Soleris / AOAC 2003.01
	Salmonella**	Negative**	BAM Ch. 5 /AOAC 2011.03
	Yeast	<1,000 CFU/g	Soleris / AOAC 997.02
	Mold	<1,000 CFU/g	Soleris / AOAC 997.02
	Coliforms	<100 CFU/g	Soleris / AOAC 991.14

^{*} FCC 9 Page 414

^{**}If Entero test results exceed 100 CFU/g then Salmonella testing (Negative/25g) is completed.